

RATIOPLANT®


SMART



THE NEW WAY OF DIRECT IMPLANTATION

Content

Introduction	03
Characteristics RatioPlant® SMART	04-05
Surgical steps SMART	06-07
Healing phase and exposure	08
Impression and prosthetical care	09
Tools and packaging	10
Implants and prosthetic parts	11
Safety, Liability and Warranty	11
Contact	12



Introduction

Our long years of experience in the field of human implantology and our expertise in the development, manufacture and testing of implants and instrument combinations provide the ultimate guarantee of extreme functionality throughout the entire HumanTech product portfolio. Given the growing need to increase the quality of human life and the dynamic market changes marked by rising pressure on costs and profit margins, factors such as cost-oriented manufacture and distribution are increasingly occupying centre stage.

HumanTech is a company group dedicated one hundred percent to the development, deployment and manufacture of implants and instruments in the medical field, and to the on-going search for ever better solutions.

From development through to the finished product and customer service – all from a single, reliable source.

The RatioPlant® implants are manufactured, packed after the up to date guidelines in our house and brought directly to our customers to the dispatch.

The variety of the RatioPlant® implant lines offers a broad range of clinical solutions, up to reconstructions of single crowns, screw connected or fixed cemented bridges and partial or full dentures. Furthermore you can use RatioPlant® of implants in all situations, from the simplest to the most difficult, surgical and bone augmentation procedures.

They are manufactured from biocompatible quality titanium and by their blasted and etched surfaces on the state of the art.

All RatioPlant® implants fulfill the highest international standards. We are certified according to DIN EN ISO 13485 as well as appendix II of the guideline 93/42 EEC.



Packaging

- All RatioPlant®-Implants are available in a new tube packaging, which is packed in a separate blister - user-friendly, safe and double sterile.
- This package ensures an easy take out with the adapter directly from the tube during surgery.
- Patients labels with all relevant data enable an easy documentation of inserted implants.



Safety, Liability and Warranty

Safety

The RatioPlant®-implant system may be used only under the guidance and recommendation of the HumanTech Germany GmbH. The use of components which are not corresponding original components to the system will impede the functionality and exclude our liability. Guidance on the use of products made verbal and in demonstration events. It corresponds to the current state of knowledge at the time of distributing our products. This does not absolve the user from his obligation to the individual product in each case before the proposed use on its suitability for the intended purpose to verify.

The processing and application of the products is up to the responsibility of each user. The liability for damage resulting from the use and application of the product is excluded. As part of our general business conditions we confirm the product quality of our products with CE certification, according to the current state of science and technology.

Dispensing

The products are delivered only to dentists, doctors, surgeons, dental technicians, dental clinics and dental laboratories.

Replacement

The withdrawal of the products can only be done in the course of an exchange. Condition for redemption of goods:

1. Two years before the end of sterility
2. Undamaged, optically modified and original packed.

RatioPlant® SMART



Description

RatioPlant® SMART implants are enossal implants available in various lengths and forms. They are inserted surgically in the bone of the upper and/or lower jaw to anchor functional and aesthetic oral rehabilitation in partially and fully anodont patients.

Prosthetic treatment may involve individual crowns, bridges and partial or full dentures, which are connected to RatioPlant® SMART implants through appropriate elements.


The RatioPlant® SMART implant system includes surgical, prosthetic and laboratory components and instruments.

In principle there are no preferred areas of use for the RatioPlant® SMART implant system.


Advantages

- **Immediate implantation** without immediate loading after tooth extraction
- **Maximum stability** more than twice as high as with other implantation systems
- **No bacteriological contamination** causing peri-implantitis by gap or pumping effect
- **Prevention of bone resorption** through stimulation of bone cells in the cancellous area due to optimised distribution of masticatory force
- **Double 3D anchoring** for better bone growth and anti-twist safeguard
- **Patented spacer technology**
- **Simplification** of typical treatment course
- **Ergonomic design** of surgical kit allows an immediate and failure-free implantation

Characteristics



Stability



RatioPlant® SMART


Surgical steps - example for immediate implantation after extraction

Surgical phase

1. Extraction of non-viable tooth and preparation of alveolus -see drilling protocol
2. Widening of alveolus with SMART final drill
3. Prepared implant bed
4. Removal of implant from sterile pack directly with implant inserter. Insertion of implant with implant inserter (A)
5. Before applying the spacer, determine the size of the spacer and the body with the measuring gauges (depth gauge SMART 085/095/105). In this case, it must be ensured that approx. 1 mm clearance between the measuring post and the antagonist is provided occlusally (B). The measuring gauges are inserted into the implant as far as the stop. If necessary, lower the implant position.
6. Insertion of spacer with spacer inserter and finger-tight screwing in (C)
7. Removal of spacer inserter
8. Optionally the bone can now be built up with suitable material to improve bone stability.
9. Before the operation wound is closed, the spacer and the bone replacement material can be covered with a suitable membrane.

Note:

When removing the sterile parts of the packaging tube, make sure that after unclip the two covers, the opening facing up to prevent falling out of parts!



Surgical steps - example for late implantation

Surgical phase


1. After exposure carry out pilot drilling with the triangle drill - see drilling protocol
2. Widening of alveolus with twist drill 2.8 and 3.5.
3. Final preparation of alveolus with SMART final drill. Optional with the countersink 5.0 at hard cortical bone.
4. Prepared implant bed. Removal of implant from sterile pack directly with implant inserter. Insertion of implant with implant inserter.
5. Before the spacer is inserted, determine the size of the spacer and the body with the measuring gauges (depth gauge SMART 085/095/105). In this case, it must be ensured that approx. 1 mm clearance between the measuring post and the antagonist is provided occlusally (B). The measuring gauges are inserted into the implant as far as the stop. If necessary, lower the implant.
6. Insertion of spacer with spacer inserter and finger-tight screwing in and removal of spacer inserter.
7. Optionally the bone can now be built up with suitable material to improve bone stability. Before closure the spacer and the bone substitute material must be covered with a suitable membrane.
8. Closing of operation wound and attend the healing time.

Note

- The following descriptions are not sufficient for the immediate application of the RatioPlant® Implant System. We recommend the briefing into the handling of the RatioPlant® Implant System by an experienced surgeon. Fundamentally the RatioPlant® Implant System only should be used by trained dentists, implantologists and dental technicians.
- Methodological errors can lead to the loss and damage to the peri-implant bone. Processing and application of the products are beyond our control and are under the responsibility of each user. Any liability for damage which caused in this case is excluded.
- Please also see our information on safety, liability and warranty on page 11.



RatioPlant® SMART



Healing phase and exposure


10. Healing of implant according to patient's situation (3-6 months)
11. Exposure by scalpel. Then unscrew the spacer with the body inserter.
12. Carefully insert the rasper instrument into the implant and activate the existing bone with light rotating movements.
13. Depending on the depth of the implant position, 3 different bodies are available with 8.5mm, 9.5mm and 10.5mm lengths of the crestal part. Before inserting the spacer, measure the size of the body using the gauge (depth gauge SMART 085/095/105). First remove the healing cap from the tube package and then remove the corresponding body with the body inserter (C) from the sterile package and screw it into the implant. Screw it into the implant and tighten with the torque ratchet and a torque of at least 40Ncm.
14. Then place the healing cap for forming and healing on the body, fix it finger-tight with the cover screw and allow to heal for a further 6-8 weeks for bone integration.

Bodies

To compensate different bone depth to implant following Bodies are available:




Length of Body in mm	8.5	9.5	10.5
Depth of implantat bed in mm	8.5	9.5	10.5
Total length of complete implant in mm	13.5	14.5	15.5
Bone heigth over implant in mm	1.9	2.9	3.9
Gingiva area in mm	1.5	1.5	1.5



Impression and prosthetic treatment

15. After integration of the body remove the cover screw, take off the healing cap and put on an impression cap. This will be in the correct position when it is perceptibly engaged on the body.
16. Then fill the impression tray with suitable impression material, coat the area around the impression cap likewise and apply the impression tray.
17. After the impression material has hardened, take off the impression. After impression taking, put the healing cap back on and fix with the cover screw.
18. Reposition the lab analogue in the impression cap and fabricate the master model. It is recommended that the mucosal section be made from soft silicon material. There is a wide range of semifinished parts available for fabricating prosthetic solutions. For further solutions please ask our technical adviser or your local dealer.
19. To prepare an individual crown for example, a plastic abutment SMART S hex(red) can be used as a base, which can be cemented on or occlusally screwed after casting and finishing. For bridges or other connected structures, like bar constructions, the plastic abutment SMART S(White) is used.
20. Before definitive fitting of the prosthetic parts, clean and drain the supragingival part of the body. Then either cement the prosthetic part as appropriate or fix it occlusally with a suitable prosthetic screw. Occlusal screwing is normally employed in the molar and premolar area and must be provided for in the metal framework.
21. Alternatively a digital impression can also be taken with the scan body or the cast can also be used to fabricate a CAD-CAM denture indirectly on the master model.



Tool	Dreikantbohrer triangle drill 2.1	Spiralbohrer twist drill 2.8	Spiralbohrer twist drill 3.5	Finalbohrer final drill	Versenker countersink optional
Durchmesser diameter Ø	2.1 mm	2.8mm	3.5mm	4.7mm	5.0mm
Drehzahl/RpM	900-1200	400-700	400-700	400-600	200-500

★ Anzuwenden bei Spätimplantation an D1 und optional bei D2 Knochen / Use at late implantation in D1 and optional in D2 type bone!

Während des Bohrvorgangs auf ausreichende Kühlung achten! / During the drilling process to ensure adequate cooling!

Tiefenmarkierungen am Finalbohrer entsprechend der Implantatbettiefe bei 8.5, 9.5 und 10.5mm / Depth markings on final drill according to the implant bed lengths of 8.5, 9.5 and 10.5mm

Um einer Schädigung des Knochengewebes vorzubeugen, ist die abgebildete Bohrfolge einzuhalten! / To prevent damage of the bone tissue, the imaged drilling sequence is observed!

MEDICAL DEVICES


HumanTech Germany GmbH - Gewerbestrasse 5 - D-71144 Steinenbronn - Germany - Telefon: +49 (0)7157-5246-71 - Fax: +49 (0)7157-5246-33 - Email: info@humantech-solutions.de

Tools

ratchet torque	5012303002	
screwdriver hex ratchet short	5012301003	
screwdriver hex ratchet long	5012301005	
screwdriver hex motor short / ISO	5012301001	
screwdriver hex motor long / ISO	5012301002	
triangle drill 21	5010315341	
twist drill 28	5010328374	
twist drill 35	5010335377	
countersink 5.0	5010335377	
final drill SMART	5003907020	
implant inserter ratchet short	5003907009	
implant inserter ratchet long	5003907010	
implant inserter motor short	5003907011	
implant inserter motor long	5003907012	
body inserter ratchet short	5003907013	
body inserter ratchet long	5003907014	
body inserter motor short	5003907015	
body inserter motor long	5003907016	
depth gauge SMART 085	5003907024	
depth gauge SMART 095	5003907025	
depth gauge SMART 105	5003907026	
rasper short	5003907019	
rasper long	5003907023	

Implants and prosthetic parts

Implant SMART 42-065	5003142065	
body SMART Ti 0 SH85 inkl. healing cap inkl. cover screw	5003114085	
body SMART Ti 0 SH95 inkl. healing cap	5003114095	
body SMART Ti 0 SH105 inkl. healing cap	5003114105	
transfer cap SMART S	5003174201	
lab analog SMART S	5003184201	
plastic abutment SMART S inkl. prosthetic screw	5003193401	
plastic abutment SMART S hex inkl. prosthetic screw	5003193411	
bar connector SMART Ti S inkl. prosthetic screw	5003193301	
scan connector SMART peek S inkl. prosthetic screw	5003203301	
prosthetic screw SMART	5003190001	



Sales Europe

HumanTech Germany GmbH

Gewerbestr. 5
D-71144 Steinenbronn

Germany

Phone: +49 (0) 7157/5246-71
Fax: +49 (0) 7157/5246-33
info@humantech-solutions.de
www.humantech-solutions.de

Sales Middle East

HumanTech Med. Sag. Tic. Ltd.

İkitelli OSB Tümsan 2. Kısım
C-Blok No: 47
TR-34306 Başakşehir İstanbul

Turkey

Phone: +90 (0) 212/485 6675
Fax: +90 (0) 212/485 6674
info@humantech.com.tr
www.humantech-solutions.de

Sales Latin America

HumanTech Mexico, S. DE R.L. DE C.V.

Rio Mixcoac No. 212-3
Acacias del Valle
Del. Benito Juárez
C.P. 03240 Mexico, D.F.
Mexico

Phone: +52 (0) 55/5534 5645
Fax: +52 (0) 55/5534 4929
info@humantech-solutions.mx
www.humantech-solutions.de

Other Countries

HumanTech Germany GmbH

Gewerbestr. 5
D-71144 Steinenbronn

Germany

Phone: +49 (0) 7157/5246-71
Fax: +49 (0) 7157/5246-33
info@humantech-solutions.de
www.humantech-solutions.de

